



# Pandemic respiratory protection masks

## Nora F CPA

Article number Nora F-009

Made in Germany/Saxony

- ✔ Corona SARS-Cov-2 pandemic respiratory protection masks (CPA) according to ZLS (Central Office of the German Federal States for Safety engineering) testing principles (without CE marking)
- ✔ tested based on EN 149:2001+A1:2009\*
- ✔ with ear bands for perfect fit on the face
- ✔ our flat fold design enables simple and space-saving storage before use
- ✔ the integrated nose clip offers a perfect fit and a secure seal





## Features

- ✔ mask made of spunlaced nonwoven
- ✔ three different functional layers
- ✔ with soft elastic ear bands
- ✔ with nose clip
- ✔ very light (5,5 g mask)
- ✔ individually packed in foil
- ✔ outer layer: stability & separation
- ✔ middle layer: ensures the filtration of the aerosols (ultra-fine particles)
- ✔ inner layer: ensures moisture transport and an optimal wearing

### Composition mask:

outer layer: 100 % PET  
middle layer: 100 % PP Meltblown  
inner layer: 100 % fine fiber (PET/PA)

### Composition ear bands:

Polyester/Spandex, very soft

### Composition nose band:

metal bracket with plastic coating

### Durability :

2 years in original packaging

### Storage:

Store at room temperature in a dry and ventilated room

PROPERTY	TEST METHOD	UNITS	TARGET	VALUES
Weight	Internal Test	g		5,5
Respiratory resistance at 95 l/min (inhale)	*Test clause of EN 149 8.9 only	mbar	≤ 2,4	≤ 2
at 160 l/min (exhale)			≤ 3	≤ 3
Aerosol Penetration Paraffin oil at 95 l/min	*Test clause of EN 149 8.11 only	%	≤ 6	≤ 3
Aerosol Penetration NaCl at 95 l/min	*Test clause of EN 149 8.11 only	%	≤ 6	≤ 3

Version 1.1 from 2020/09/09

Norafin Technologies GmbH

Gewerbegebiet Nord 8

09456 Mildenau | Germany

Tel. +49 3733 5507 0

info@norafin.com

[www.norafin.com](http://www.norafin.com)

